

Velashape II body slimming machine

User manual

装
订
线



Thank you for buying our products.
Please read this instruction carefully before use, and operate in strict accordance with the instructions.
After reading, please keep the instruction for future reference.

CONTENTS

Chapter 1 Introduction to the instrument	1
1.1 Technical principles	2
1.2 Instrument parameters	2
1.3 Product List.....	4
Chapter 2 Instrument Operation	5
2.1 40K Cavitation handle operation	5
2.2 Large Vacuum Roller operation.....	5
2.3 Medium vacuum roller operation.....	5
2.4 Vacuum face lifting handle operation.....	5
Chapter 3 Interface system Introduce	6
3.1 40K operation interface	6
3.2 Large Vacuum Roller operation interface	6
3.3 Medium vacuum roller Operation Interface.....	7
3.4 Vacuum face lifting handle operation interface.....	7
Chapter 4 Operation Process and Precautions	8
4.1 Body operation steps	8
4.2 Facial operation steps	8
4.3 Operating precautions	8
4.4 Operational contraindications.....	9
4.5 Routine maintenance and repair of the instrument	9

Chapter 1 Introduction to the instrument

The vertical Velashpe 2 body slimming and shaping machine is a multifunctional operating system integrated machine that combines the body and the face, which can improve problems such as body obesity and facial aging. The instrument acts on the human subcutaneous fat cells through safe high-frequency energy waves. The same frequency resonance is generated inside fat cells to generate powerful shock waves, free radicals and high-speed microjets. Safely decompose the fat cell wall from inside to generate cell debris and free fat. Cell debris is engulfed by phagocytes, and free fat is metabolized to the liver via lymph and blood. It can effectively improve the curve of the legs, enhance the sagging buttocks, reduce the abdomen, enhance the firmness of the skin, and improve the stubborn hardened cellulite. Compared with traditional treatment methods, the advantage is that it is non-invasive and does not cause any damage to the body. It is a better non-surgical shaping method; the negative-pressure radio frequency head of the face can stimulate collagen hyperplasia and elastic fiber reorganization, effectively improving facial relaxation and Wrinkles, radiofrequency heating of subcutaneous fat, eliminate double chin obesity position, plus vacuum negative pressure to dredge the lymph, immediately improve tightening and anti-aging effect.



- 1 40KHZ cavitation head
- 2 Large Vacuum Roller head
- 3 Screen
- 4 Medium vacuum roller head
- 5 Vacuum rf face head
- 6 ON/OFF

1.1 Technical principle



(1) 40K cavitation principle

The cavitation with a frequency of 40000HZ is emitted by an ultrasound that gathers ultrasound, and the energy wave enters the subcutaneous fat cells of the human body, and generates the same frequency resonance inside the fat to produce powerful shock waves, free radicals, and high-speed microjets. Decompose to produce cell debris and free fat. Cell debris is engulfed by phagocytes, and free fat is metabolized to the liver via lymph and blood.



(2) large roller Principle

Using vacuum suction and a unique roller design, the skin and subcutaneous muscle tissue are lifted up, and then the wave deep massage is performed on the body through the role of the vacuum suction front and rear rollers to suck, pinch and roll the skin and fat tissue. Wait for mechanical aerobic exercise. The tri-color lamp functions red light to accelerate blood circulation, blue light to calm the skin, and green light to promote absorption. Deep massage accelerates the metabolism of triglycerides and fatty acids decomposed by 40K strong sonic waves, and is excreted through the lymphatic circulation. This vacuum deep massage can increase the lymphatic detoxification function by four times and the skin blood circulation by four times, making the body more effective. Consume fat accumulation to achieve slimming effect.



(3) Middle Roller Principle

Using vacuum suction and a unique roller design, the skin and subcutaneous muscle tissue are lifted upwards, and then the wave deep massage is performed on the body by the role of the vacuum suction front and rear rollers to the skin and fat tissue perform mechanical aerobic exercises such as suction, pinch and roll. Infrared light accelerates blood circulation. Deep massage accelerates the metabolism of triglycerides and fatty acids decomposed by 40K cavitation, and is excreted through the lymphatic circulation. This vacuum deep massage can increase the lymphatic detoxification function by four times and the skin blood circulation by four times, making the body more effective. Consume fat accumulation to achieve slimming effect.



(4) Facial vacuum handle

Use vacuum, radio frequency and infrared to improve aging problems such as facial sagging and more collagen loss, strengthen blood circulation to the face, activate sleeping cells, promote collagen regeneration, solve excess fat on the face and relax sagging and wrinkles, help the face lift. Firms, lightens wrinkles, and lifts contours.

1.2 Technical parameter

① Machine parameters

Name	V9 Velashape II body slimming machine
Input voltage	220V/50HZ 110V/60HZ
Size	120cm*142cm*66cm
Power	800W
Fuse	220V-10A/110V-20A

② 40K cavitation

Name	40k cavitation handle
Operation voltage	36V
40K output frequency	40KHZ
Power	54W
Effective area	33cm ²
Energy density	1.62W/cm ²
Energy output	pulse
Output Level	8

③ Facial vacuum handle

Name	Vacuum face lifting handle
Vacuum pump working voltage	220V/110V
RF working voltage	DC36V
RF output frequency	5MHZ
Out put level	10

④ Middle vacuum roller handle

Name	Middle vacuum roller handle
Vacuum pump working voltage	220V/110V
Motors working voltage	DC5V
RF working voltage	DC36V
RF Frequency	5MHZ
RF output level	10
Infrared wavelength	650nm

⑤ Large Roller Handle

Name	Large vacuum roller handle
Vacuum pump working voltage	220V/110V
Motors working voltage	DC12V
RF frequency	5MHZ
RF output Level	10档位
RF power	60W
Vacuum pump power	270W
Infrared wavelength	650nm
Green light wavelength	560nm
Blue light wavelength	490nm

1.3 Accessories list



Host



Large Roller



Middle Roller



Vacuum Facial



40K cavitation



Power cable

Chapter 2

Instrument Operation

2.1 40K cavitation handle operation



Apply essential oils or slimming cream products to the operation site, and use 40K handle to act on the operation site in a circular motion. The operation time is set to 20-30 minutes.

Note:

- (1). During the operation of 40k cavitation, the client will have tinnitus phenomenon, which is normal, because when 40K ultrasound waves act on the skin, the sound waves scattered on the fat will be transmitted to the ears along the skin.
- (2). The operation must be close to the skin during the cavitation operation .
- (3). 40k cavitation operation can be used with essential oil or slimming cream.
- (4). 40K cavitation cannot be operated on bones, and long-term operation of ultrasound waves on bones will cause osteoporosis.

2.2 large roller handle



The large roller contains the functions of vacuum, radio frequency, four-direction roller and three-color light. If you are using a slimming cream product when operating 40K, you need to clean the residual cream in the operating part first. If it is not clean, it will cause internal air tube blockage , You need to use essential oil as a medium, and the roller needs to be completely close to the skin during the operation, avoiding the navel. Operating time is 25 to 35 minutes.

Note:

- (1). It is necessary to use essential oil as the medium when operating the roller head. It cannot be operated with milk, cream or paste products, which will cause the internal trachea to be blocked.
- (2). Due to the different textures of essential oils, it is normal for some essential oils to turn black during operation. Chemical phenomena caused by the exposure of essential oils to RF heating
- (3). The large roller is suitable for the parts with large operation surface, such as abdomen and thighs.

2.3 Medium roller handle operation



The medium roller contains vacuum, radio frequency, four-direction roller and infrared light function. If you are using a slimming cream product when operating 40K, you need to clean it with clean water first. If it is not clean, it will cause the internal trachea to be blocked. Essential oil needs to be used as the medium, the rollers need to be completely close to the skin during the operation, and the operation time is 25 to 35 minutes.

Note:

- (1). It is necessary to use essential oil as the medium when operating the roller head. It cannot be operated with milk, cream or paste products, which will cause the internal trachea to be blocked.
- (2). Due to the different textures of essential oils, it is normal for some essential oils to turn black during operation. Chemical phenomena caused by the exposure of essential oils to RF heating
- (3). The medium-sized roller is suitable for small parts such as calves and arms

2.4 Facia vacuum handle



The facial vacuum handle of the face contains vacuum, radio frequency, and infrared rays. Essential oils or lubricating, volatile massage creams are required for operation. Paste products are not available. The operation needs to be completely close to the skin.

Note:


Because the facial skin is relatively sensitive, the vacuum knob should be adjusted to the minimum position before operation. The point on the knob on the back of the fuselage is the minimum value of 20 scale, and slowly increase according to the strength of the skin. In order to avoid excessive negative pressure, direct operation will cause skin discomfort.

Chapter 3 Interface system Introduce




Connect the power and turn on the key,
Press **ENTER** into the operation interface

3.1 40K operation interface



Blasting fat 40KHZ shows as the energy intensity, press **Blasting fat 40KHZ** after it turn red, You can use the up and down arrows to increase or decrease the gear, the energy is up to 8 gear, 40K operation head time setting interface, the working time is 30 minutes by default, press **Time** when it turn to red, use the arrows to adjust the operation times.

3.2 Large vacuum Roller Operation




The direction of the 4 rollers can be switched arbitrarily according to the user's feeling **Time**. The default working time is 30 minutes, press **Time**, when it turns into red to increase or decrease the time by pressing the arrow. **R.F** the default rf energy is level 1, and 10 is the highest. Press the right and left arrow to control the output energy. Press **R.F** to turn the rf on. There are 3 air pressure output, mode A: suction 1.5 seconds and release 1.5 seconds, mode B: suction 1 second and release 1 second. Mode C: suction 0.5 second and release 0.5 second.

Press **Laser OFF** to turn on or turn off the air pressure, adjust **Laser RED** by pressing the up and down arrow to increase or decrease the air pressure output according to user's tolerance.


Press **Laser RED** to change the light output, red, green, and blue light, or **Laser RGB** 3 lights could be work together. Turn light on or off by pressing **Laser OFF**.

3.3 Middle Roller handle



press **R.F (w)** choose the middle roller operation control panel, time adjust **Time**, RF adjust **R.F (w)**, turn on the RF, choose the mode, turn on the Vacuum, press the button on the head of the handle to turn on the roller, 4 lights is perform as 4 roller working modes, press one time to switch the modes, long press it 3 seconds to turn on the infrared light, and press vacuum **Vacuum** and rf button **R.F** on the main screen to start working. when the red light turn to the last mode, press one second to stop the roller working, press the button 3 second to turn the infrared light.

3.4 Vacuum face lifting handle



press **R.F (w)** to choose the operation control panel, time adjust **Time**, RF energy adjust **R.F (w)**, turn on the rf and choose the mode, vacuum and infrared light too, **R.F** like this when it turn red you can start to work. If the client want the vacuum stronger, you can use this button behind the machine to adjust the vacuum.

Chapter 4 Operation Process and Precautions

4.1 Body operation steps

- (1). Clean the operation area with warm water, apply the medium essential oil or slimming cream products evenly in a circular motion, and turn on the 40K cavitation. It is recommended to start the operation from energy 3, adjust it according to the energy intensity of the customer, and adhere to the skin Operate for 20-30 minutes.
- (2). Clean the residual paste during the 40K cavitation operation, apply the slimming oil evenly in a circular and soothing technique, and turn on the large roller. It is recommended that the RF energy intensity starts from 3 and starts to operate against the skin for 20-30 minutes.
- (3). After the operation of the instrument is completed, use the soothing and detoxification technique to drain the hands alternately in the direction of the inguinal lymph, massage the operation part for 10 minutes, clean with warm water, and the operation is completed.

Treatment recommendations: 12 times a treatment course, at the first 3 times, one time one day, the next 9 times, every 2 days for one time.

4.2 Face operation steps

- (1). Cleansing products clean the face
- (2). Facial anti-aging enhances firming oils, evenly apply to the face, turn on the handle and set up the parameters, it is recommended that the RF energy intensity starts from 3, adjust according to the energy intensity of the skin to bear, start the operation with the lifting method, the operation time the whole face 15-20 minutes.
- (3). applying mask after the operation, cleaning the handle, finished

Treatment recommendations: 10 treatments per session, once a week.

4.3 Operation notes: (Please clean the operation head with warm water before use)

- (1). Use a vacuum handle and roller do not use gel or cream products. Essential oils are preferred.
- (2). When using the vacuum head handle and roller handle, it needs to be close to the skin and wait for the skin to move.
- (3). 40K cavitation needs to be close to the skin when used.
- (4). When the operation head is not in use, try to close the related operation buttons.
- (5). After using the device for a period of time (about 1-2 hours), it is best to shut down and stop working for 10 minutes, and then restart the instrument.
- (6). Vacuum handle and roller handle. After use, massage hands can be used to alternately wipe and drain to lymph nodes near the body to facilitate better metabolism of body toxins.
- (7). For slack fat, the use time of the roller handle can be extended by 5 minutes: for stubborn fat, the use time of the 40K cavitation handle can be extended by about 5 minutes.
- (8). Customers can be advised to drink a cup of warm water before the operation, which is conducive to speeding up the blood circulation of the skin and the discharge of toxins from the body.
- (9). After the operation, carbohydrate foods (such as rice, noodles, desserts, fruits, etc.) and fats should not be ingested as much as possible during this period, so that the body burns the decomposed fat as energy.
- (10). Edible tomatoes, cucumbers, boiled vegetables, eggs and other foods

4.4 Operation contraindications

- (1). It is recommended to avoid skin wounds or scars;
- (2). Use hormone drugs with caution;
- (3). Physiological period of pregnant women and women;
- (4). Inside the body to metal sheets, etc. (such as pacemakers, etc.)
- (5). Those with extremely weak bodies
- (6). Within three months after a woman wears a birth control ring
- (8). Women should not operate the abdomen within three months after delivery, and the abdomen should not be operated within six months after caesarean section
- (9). Female lactation period
- (10). Skin during allergies

4.5 Daily maintenance and repair of the instrument

- (1). All handles of the device are made of ABS resin material, and the instrument should not be wiped with alcohol everyday.
- (2). The device is a liquid crystal touch screen. When tapping, try to switch with your fingers, and do not touch with sharp objects.
- (3). When using the accessory, try to handle it lightly, and do not drop it by gravity: try not to pull the output line during operation.
- (4). When each accessory is not in use, try to shut down the work as much as possible: after the instrument is operated for 1-2 hours, try to pause the instrument for 10 minutes and then restart the work.
- (5). Please clean the instrument with warm water every time you use it. Wipe the entire handle and handle operation surface with a wet wipe. Each handle must be wiped clean. There should be no residues such as product media, and the power should be properly stored after disconnecting the power supply.
- (6). Before using the instrument, you must use a plug with a grounding foot, and make sure that the power socket of the instrument is truly grounded.
- (8). In order to ensure the therapeutic effect of the instrument and ensure the normal service life, please use the specified accessories provided or recommended by the original manufacturer.
- (9). The instrument should not be placed in a humid place or near a water source, and the instrument should not be directly exposed to sunlight.
- (10). Do not put the instrument close to a strong heat source, because it may affect the life of the instrument and normal use.
- (11). Use the current loop arm sleeve in winter or dry weather. To avoid poor contact, please keep a certain degree of moisture in the place where you touch the skin. Wipe it with water-based products or water to achieve the desired effect.
- (12). Water, oil and other media must not penetrate the inside of the instrument, and can not be struck, knocked or dropped to avoid damage to the instrument.
- (13). The warranty period is one year, under the premise of normal use of the machine during the warranty period, it is responsible for free maintenance and replacement of parts, except for consumables (circuit board, air pump, motor, display)
- (14). All accessories are prohibited from violent collision (especially 40K cavitation handle) and drop.
- (15). Provide lifetime maintenance and technical support.