

SPA100 MULTIFUNCTION SKIN SPA SYSTEM



IMPORTANT: before starting the equipment correctly, please follow the installation instructions. Any damages caused by incorrect operations will void the warranty.

SAFETY INFORMATION

Thank you for choosing our products! This user instruction is specially

designed for the safety of you and others. Please read the instructions carefully before installing or using the machine and follow them when using.

1.1 Electrical Safety

-Before installing the machine, make use that your local electric network complies with the voltage and frequency indicated on the tag on the rear panel. Otherwise, the machine may be damaged in case of excessive current or voltage rise.

-Never touch the machine with damp or wet hands or feet.

-It is dangerous to change or try to modify the function of this machine in any way.

-Do not attempt to operate this unit if the power cord has been frayed or broken

-Do not spill water or other liquids into or on to your unit.

-During long periods of non-use, disconnect the unit' s main power.

1.2 Safety Rules When Installing

-Place the machine on a stable and secure base. It should not be subjected to any jolts or vibrations.

-Never expose the unit to too high ($>40^{\circ}\text{C}$) nor too low temperature ($<5^{\circ}\text{C}$).

-NEVER EXPOSE THE UNIT TO RAIN OR ANY OTHER HUMIDITY.

-Do not expose the unit to direct sunshine or install the unit close to a source of heat (central heating, radiator...); otherwise the screen may be damaged.

-Do not cover or block the ventilator behind the unit.

1.3 Safety Rules When in Use.

-Well clean areas to treat before treatment.

-Never touch the machine with damp or wet hands or feet.

-Never place the cavitation handpiece on the liver, kidney area.

INTRODUCE OF THE FUNCTION

1. Hydro Dermabrasion

Hydradermabrasion (Hydra Facial) – the latest in skincare technology. Hydradermabrasion utilizes the natural healing powers of water and oxygen to effortlessly exfoliate the skin without the use of hard crystals or abrasive textured wands, generating deeply hydrated and healthier looking skin.

Application

1 Rejuvenate sun damaged skin-face, neck, shoulders, back,


arms and legs. 2 Reduce age spots

3 Minimize of blotchy skin coloring

4 Reduce acne and superficial scars from past injury 5 Remove blackheads and white heads

6 Reduce oily skin

7 Improve overall skin health



2. Diamond Dermabrasion

In this type of treatment, dead skin cells are exfoliated and vacuumed off the face with a microdermabrasion machine, that uses microdermabrasion tips covered in tiny diamonds. This procedure is non-surgical, known to be most effective on deep acne scars or old stretch marks.

Besides, the diamond tips together with the strong suction make it possible to exfoliate the skin more thoroughly than other

procedures. With the flexibility of using different suction levels, it also increases blood flow to the skin, and that in turn helps to produce collagen.

Features

- 1.Suitable for all skin types
2. Provides non-invasive, non-surgical procedure options
- 3.reliable for skin issues like wrinkles, fine lines, acne scars, uneven skin tone, and stretch marks

3. Spray treatment Mechanism of spray oxygen

Operating steps

- 1.Clean skin deeply.
- 2.Click the “Spray gun ” on as Fig.2 shows,
- 3.You can set the operation time .
- 4.Mix 2 ~ 3 drop SPA essential oil diluted with purified water into small cup of the sprayer,water should be about 2 / 3 of the cup, and then tighten the lid tightly. Oxygen sprayer can only with liquid products, other states product will cause pinhole jam.
- 4.Beautician hold the spray gun, press the control button to release oxygen, And pull backwards to spray water mist.
- 5.Spray oxygen from top to bottom starting from the forehead to the neck, careful not to spray into the eyes or nose, in the spray

oxygen process, the beautician can appease, lift up, stretch wrinkles with the other hand, more helpful for nutrition products penetrate the skin, strengthen skin elasticity, spray rate should depend on the guests' absorb state of oxygen.

Cautions

The product should be liquid when using oxygen sprayer.

2. Avoid spraying into mouth or nose when spraying oxygen.

3. Clean the small cup and the pipe inside the sprayer every time use it, to avoid jam.

Clean method:

pour the purified water into the small cup, and press the control button and then pull back towards.

1. Supply nutrients and moisture for the skin rapidly, smooth skin, lightening stains, improve chlorosis, dull conditions, make skin rosy, white, transparent, tender, tightening indeed, with special spa oils will be better.

4. BIO Microcurrent

Microcurrent is a low level of electricity that mirrors the body's own natural electrical impulses. It is a non-surgical,

non-invasive cosmetic technique that safely and effectively improves the health and appearance of the skin. When used in conjunction with specialized products and manual manipulations, these tiny microcurrent impulses encourage your body's currents at a cellular level. The signs of aging are greatly reduced, while skin tone and elasticity are dramatically improved.



Operation Instruction

Make sure the power line has been installed correctly; make sure the hand piece is well connected with the equipment; Ensure the circuit breaker (if any) at the back of the machine is pressed to "on" . Please find the figure 1



Figure. 1

1. Select Hydro Dermabrasion/Diamond Dermabrasion function,

(find the figure.2), then adjust the “” and “” for the water level and vacuum, working time setting

  , start working.




Figure.2

2. Select Bipolar function(find Figure.2), then adjust the



for the intensity, working time setting, start working.

3. Select Spray gun function(find Figure. 2), then adjust



for the working time, start working.

Specification



Hydro Dermabrasion



Diamond Dermabrasion



Spray gun



Bipolar

Tips:



Model	SPA100
Input Voltage	AC100V~240V, 47~63Hz
Output Voltage	150VA
Max output vacuum pressure	100Kpa
Handle	Hydro Dermabrasion/Diamond Dermabrasion/Spray gun/Bipolar
Machine Power	250W of German pump
Packing Size	62*48*42cm
NW/GW	24/25KGS
MOQ	1PC