



PORTABLE EYES BEAUTY MACHINE (DEVICE)

Instructions for use

Catalog

Product overview	01	Precautions	07
Instrument parameter	02	Instrument maintenance	08
System interface	03	Packing list	09
Operation procedure	04	Warranty card	10–11
Operation method	05		
Using taboo	06		

Product overview



Welcome to use BB-eyes Eye beauty machine,It' s a kind of new version eye beauty machine. With the technology of intelligent organic eye system (CLIP) and biological current (EMS) to gain beautiful eyes. Please read this manual carefully before use and use the beauty equipment in strict accordance with the guidelines in this manual.

Instrument parameter





Product Name: BB*eyes	Model Number: Z-1805
Voltage:100~240V	Power: 2.8W
Input voltage: 9V	Operating probe: 1
Microcurrent power: 1.3W	Demensions: 18*4.2*4cm
Executive standard: GB4706.1-2005	Feature: EMS/CLIP
Package weight: 2KG	widths and balance : 14US
Package size: 31.4*20.07*7cm	

System interface



Operation procedure

-  Clean eye around before using the machine (If you have make up , you need to use make up remover to remove the make up then use facial cleanser to clean)
-  Connect the power cable, There is a round button above the power supply. It is the button to open the intelligent organic eye system (CLIP) and biological current (EMS), press these 2 buttons and hold for 2 seconds or more separately to turn off the machine.
- CLIP has 3 frequency adjustments, L is low speed, M is medium speed, H is high speed, It shows L low speed when open the machine, You can press the round button to adjust the speed. Biological current (EMS) share the same adjustment, The light will always blink when the micro current at the time of output. ;
- During the treatment , put the eye essence on the treatment area, it is recommended to add essence immediately when essence has been absorbed. It is recommended to apply a little bit of force to the skin during operation. Operate in a uniform lifting speed.
- Using wet wipe to wipe off the remaining essence on the probe after using the instrument ;
- No need to clean eye around after treatment, Can apply eye mask or not after treatment (It is not recommended to clean eyes with hot water , warm or cold water can do.)

Operation method



EYES



Protein ball

(1) Probe lifting from inner eye corner to eye around, probe lifting with CLIP speed from Intraocular angle to Extraocular angle, from extraocular angle lifting to temple, from temple lifting to hairline. It is recommended to avoid the upper eyelid position (upper eyelid close to eyeball, not comfortable)

Used with special-purpose products can get better effect. Special purpose products for eyes : Anti-age eye essence+Anti-age eye protein balls

Attention:

About eye treatment usually 10–20 minutes it is not recommended for more than 20 minutes, Normal skin is recommended to use it once a day for the first 5 times. After the first 5 times , then can do once every 2 days. People with sensitive eyes is recommended to operate once every 2–3 days.

Eye: Before operating , open the eye protein balls , Squeeze 1ML eye essence from protein ball to dissolve protein ball, After dissolving, the essence is translucent and slightly viscous. smear the essence on the eye below, then start treatment.

Noted: When feel micro hemp during treatment , that is the CLIP is lift and squeeze eye around , It is normal for eyes to turned red after operation.

Using taboo

People with the following symptoms should be caution , please consult doctor or professional before using .

- ❶ Eyes in allergy period is not recommended to use , it is recommended to use after one week of the recovery of allergy.
- ❷ Eye around has been injected or filled Hyaluronic acid is recommended to use after 30 days.
- ❸ Women pregnant is recommended to open only CLIP low speed to use , Human who are sensitive is not recommended to use.
- ❹ Eye around with wound or pustules is recommended to use after the recovery.
- ❺ Patients with severe skin diseases or transfections should be used with caution

Precautions

- ❶ It is recommended to apply water mask for 2–3 times a week.
- ❷ Less smoke, less drink alcohol, less stay up , drink more water.
- ❸ Keep the operating probe in a wet state during operation, in the same time, the probe should close to skin.

Instrument maintenance

- 1** Clean with water every time after use the instrument and save it properly.
- 2** To assure the voltage is adapted, If the local voltage is not stable , We recommend users to increase the power matching power supply. ◦
- 3** To assure the warranty. Please use the original spare parts provided by original factory.
- 4** The instrument is not recommended to be placed in a wet place or near a water source. nor put exposed to the sun. ◦
- 5** Do not bring the instrument close to a strong heat source, in case of influence the warranty and normal use of the machine
- 6** The warranty is one year, during the warranty period , if the machine is broken due to no human damage. Factory bear the fees of repairing and changing the spare parts. (clients need to bear the shipping fee)
- 7** When warranty expired , factory can help repair the machine (only charge cost of the spare parts). ◦

Packing list



Name	Cable: one
Instruction sheet: ONE	Certificate: one

Warranty card

Warranty Card	
To provide the better service , please kindly fill in the user information :	
Customer name:	Purchasing ID:
Purchasing date:	Purchasing shop:
Phone:	Wechat:
Reason:	
Warranty:	
.About all the BB–eyes products warranty shall be implemented from the date of purchase.	
.All BB–eyes products can only be repaired in location identified by BB–eyes	
during the warranty, It is suggested to show warranty card , purchasing time and certificate	

Warranty card

Noted:
<ul style="list-style-type: none">• This warranty card will not be reissued after it is lost. Please keep it;
<ul style="list-style-type: none">• This warranty card is valid only in the People's Republic of China;
From the date of purchase during the warranty period, the following situation cannot
be repaired freely , it need charge repairing fee
<ul style="list-style-type: none">• After one year warranty from the date of purchasing.
<ul style="list-style-type: none">• No warranty card , cannot prove the purchasing date and product.
<ul style="list-style-type: none">• Broken due to force majeure such as : fire, Earthquake, typhoon flood and war. Etc
<ul style="list-style-type: none">• Appearance damage such as scratches and stains generated during use
<ul style="list-style-type: none">• Broken due to personal disassemble mis repair
<ul style="list-style-type: none">• other consumable parts
<ul style="list-style-type: none">• Broken caused by general use outside the home (such as in business activities).