# Facial blackhead suction beauty machine

# → 4 suction heads → Switch → USB charging

Product accessories

# Instructions:

1, Dry your hand before using the machinein case that any liquid comes into the machine body from the vent hole on the back, which will damage the machine.

2, Make sure that the filter sponge is properly installed in the suction hole before using, and the sponge and the area to be treated are dry. be treated are dry.

3, Do not wash the machine with water.

# Operation methods Smoothingand tightening skin

1, Change to the microcrystalline head or big round holehead. Make your face clean and dry.
2, Move the machine upwards from the bottom of your chin.
3, Move the machine down the bottom of the Tzone.
4 Place the instrument on the cheek and move from the inside to the outside 5Place the instrument in the forehead area, moving from the middle to both sides

# Cleansing pore deeply

1, Change to the big round or small round head (choose according to your skin condition).
2.Aim the roundholeat the location of the blackheadsand acne to remove them. The levels can be adjusted according to your need. Method: When removing blackheads, move slowlyon the loc of the blackheads for several times. Do not move back and forth, but move in one direction. For acne, suck and then quickly unplug

#### Fading fine lines

1, Change to the oval hole head.
2, Hold the machine by hand and cleanse eye corners, around nose andmouth etc. Levels ofthe machine can be adjusted according to your need. Method: Suck the skin to be cleansed, quickly unplug, and then suck again, unplug again. Repeatfor several times, which will take effect on the skin deeply, stimulate skin'sblood microcirculation, increase skin elasticity, and prevent from growing deep wrinkles.

Installation notes of replaceable suction head and sponges







# 1, Put the filter sponge into the suction hole. 2, Press the filter cover into the suction hole. 3, Press the suction head into the nozzle.

1, The machine has lifting effect, andit is not recommended to be used on wounded, scarred skin;
2, Microcrystallinehead has microdermabrasion effect. For dry skin, use once a week; for oilymixed skin, use twice a week. Use for no more than 5 minutes each time depending on the skin condition.
3, Do not force the head on the skin, as the machine has lifting effect.
4, After using, you may feel a little irritated or the skin will be a little red, which is a normal phenomenon, andthis phenomenonwill disappearafter 10 minutes. to minutes.

5. Afterusing, make sure to moisturize your skin, prevent it from radiation, minimize skin pores and so on.

6. Before using this machine, do not use scrub cream or other exfoliating products, but steam your face over hot water to open the pores.

#### In order to use the facial blackhead suction beauty machine correctly and effectively, please make sure to read this manual, instructions and notes, and keep the manual well.

Manual

the manual well.

Product description:

Facial blackhead suction beauty machine is a multi-function beauty machine, which integratesdermabrasion, skintightening, pores, blackhead and acne cleansingfunctionsin one machine. It uses vacuum suction device equipped with microcrystalline particles torub off aging top skin layer and remove dirt inside pores, so that the pores becomefine immediately, skin becomesbright, tender, uniform and smooth. This is a non-invasive and non-irritating technology, which controls the degree of microdermabrasionbychangingthe fineness degree of the microcrystallinesuction head. With 4 different shapes of the suction heads, this product integrates 4 different effects, that is, microdermabrasion, skin tightening, pores cleansing and acne removingin it.

#### Working principle Pneumatic vacuum face lifting V shape technology:

The vacuum suction system gentlylifts and massages the skin, making the dermal tissue get enough nutrition, and the skin becomingsilkier; enhancing the permeability of the skin, so that beautyserum can penetrate into the deeper skin tissue, and thus improving skin moisture, which makes the skin becomebrighter and shinier; enhancing the skin's defense capabilities, to avoid the dirt damaging the skin, and to maintain skin firmness and elasticity, and thereby preventing skin from wrinkle growing and aging, promoting the regeneration of epidermal cells growing andaging; promoting the regeneration of epidermal cells, enhancing the skin immunity and the ability of resisting ultraviolet radiation, avoiding melanin, keeping the skin bright, and making the skin healthier; improving skin microcirculation, and therebyfading melanin spots.

# Microcrystalline probe

Product parameters:

Microcrystalline probe equipped with natural mineral microcrystalline particles, can gently exfoliate, and make the skin become more smooth and renewed. It will gently rub off rough skin surface, debris andsuck out the dirt from the skin. While removing dead skin cells, it will benefit natural renewal of cells and keep the skin smooth.

Product Name: Facial Cleaning Beauty Instrument
Battery capacity:720mAh
Voltage: 3.7V
Rated current: 230MA
Rated power: 0.8W
Charging method: USB charging
Product size: 175 \* 46 \* 47mm
Weight: 150g Weight: 150g Product accessories: 1, main machine x1 2, microcrystalline headx1, big round headx1, small round headx1, 3, filter spongex 1 package 4, USB charging line x 1 5, product manual x 1

#### Product effects



1, microcrystalline head(exfoliating): its microcrystal line head equipped with microcrystalline particles is of moderate suction intensity. With no harm to the skin , it massages the skin better, and at the same time has exfoliatingeffect, achieving skin renewal.



2, big round head (sucking out blackhead, skin tighten ing): big round head has strongsuction, and is used for sucking out blackhead andgrease. Appreciate for theskin lifting and massage effect throughsuction function of the machine, it has deep effect on underlying skin, stimulates skin's blood microcirculation, which leads to skin tightening.



3, small round head (deep cleansing): small round head has weak suction. Aim the hole at the blackhead or acne location to suck out dirt, blackheads, grease and acne. It is suitable for thin and sensitive skin.



4, oval headfading fine lines and increasing the 4, oval neadrading fine lines and increasing the elasticity: its oval design is special for lifting and massaging places, where fine lines often appear and is not easy to be cleansed, such asthe eye corners, nose sides, mouth cornersetc., so as to promote skin's blood microcirculation, increase skin elasticity and

# Areas applicable

Microcrystalline head:
any part of the faceexcept the eye corners, also rough and wrinkled
skin, such as the body / elbow and knee etc..
Big round head:
any part of the face except the eye corners.
Small round head:
sensitive skin such as eyes, also other parts of the face.
Oval head: eye corners, nose sides, mouth corners, also neck.

# Product features:

1, Interchangeable suction heads: Suction heads are interchangeableaccording to different effects as required..
2, adjustable intensity: High, medium and lowlevels are optional according to your need.
3, easy to clean: Four kinds of suction heads can be easily removed and cleaned, and with filter sponges, dirt and aging keratin are effectively prevented cominginto the machine. Sponges can be reused after cleaning.

after cleaning.
4, physical skin care: The machine is safe / convenient to operate with 4, physical skin care: The machine is safe / convenient to operatewith no edge effects or side effects. Youcan work or join in community activities immediately after use. Occasionally there will be skin redness after use, which is a normal phenomenon, andthe phenomenon will automatically disappear as soon as you stop use it. 5, small, easy to carry:Rechargeable, small and easy to carry. 6, high-capacity and rechargeable: With built-inhigh-capacity rechargeable battery of 850aMh, 5V USB portable design, the machine can be chargedwith Charge Pal, mobile phone charger of universal charging plug, or devices with USB connectionsuch as computers and portable lightsso as to achieve beauty care at home, office and outdoors.

portable lightsso as to achieve beauty care at nome, office and outdoors.

7, charging LED: When the machine is being charged, the blue light flashes in cycles. Whenthe machine is fully charged, the blue light keeps being on.

8, low voltage indicator: When the battery voltage is lower than the required voltage, the weak suction I or automatic shutdown indicates that the machine needs to be charged.

9,note: Keeppressing the switch for 3 seconds to turn off the machine no matter it is in what level.