



科瑞仪器
Kori instrument

Xinyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xinyang,Zhengzhou City,450100,China

Tel:0086-371-86573886


Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Jacketed glass reactors





科瑞仪器
Kori instrument

Xingyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xingyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Technical data 1L 2L 3L

Model	SF-1L	SF-2L	SF-3L
Size	1L	2L	3L
Motor	40W	90W	90W
Speed	60-1000RPM	60-1000RPM	60-1000RPM
Stirrer	Stainless steel with PTFE Coating 7mm	Stainless steel with PTFE Coating 7mm	Stainless steel with PTFE Coating 7mm
Maximum Temperature	200℃	200℃	200℃
Minimum Temperature	-120℃	-120℃	-120℃
Maximum Vacuum	-0.098Mpa	-0.098Mpa	-0.098Mpa
Voltage	220V	220V	220V
Package	Craft Box	Craft Box	Craft Box

Note: Jacketed glass reactor can be for heating and cooling ,but you need purchase heater or cooler additionally



科瑞仪器

>>> Kori instrument

Xingyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xingyang, Zhengzhou City, 450100, China

Tel: 0086-371-86573886

Fax: 0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130



Name	5L Jacketed glass reactor
Model	SF-5L
Size	5L
Motor	90W
Speed	60-1000RPM
Stirrer	Stainless steel with PTFE Coating 7mm
Maximum Temperature	200°C
Minimum Temperature	-120°C
Maximum Vacuum	-0.098Mpa
Voltage	220V
Package	Craft Box

Note: Jacketed glass reactor can be for heating and cooling ,but you need purchase heater or cooler additionally



科瑞仪器

>>> Kori instrument

Xingyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xingyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130





科瑞仪器

>>> Kori instrument

Xinyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xinyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130





科瑞仪器
>>> Kori instrument

Xinyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xinyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Technical data 10L 20L 30L 50L

Model	SF-10L	SF-20L	SF-30L	SF-50L
Size	10L	20L	30L	50L
Motor	90W 3:1 Gear Motor	90W 3:1 Gear Motor	90W 3:1 Gear Motor	120W 3:1 Gear Motor
Speed	60-600RPM	60-600RPM	60-600RPM	60-600RPM
Stirrer	Stainless steel with PTFE Coating 12mm	Stainless steel with PTFE Coating12mm	Stainless steel with PTFE Coating 12mm	Stainless steel with PTFE Coating 12mm
Maximum Temperature	200℃	200℃	200℃	200℃
Minimum Temperature	-120℃	-120℃	-120℃	-120℃
Maximum Vacuum	-0.098Mpa	-0.098Mpa	-0.098Mpa	-0.098Mpa
Voltage	220V	220V	220V	220V
Package	Craft Box	Craft Box	Craft Box	Craft Box

Note: Jacketed glass reactor can be for heating and cooling ,but you need purchase heater or cooler additionally



科瑞仪器

>>> Kori instrument

Xingyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xingyang, Zhengzhou City, 450100, China

Tel: 0086-371-86573886

Fax: 0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130





科瑞仪器
Kori instrument

Xingyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xingyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Technical data 100L 150L 200L

Model	SF-100L	SF-150L	SF-200L
Size	100L	150L	200L
Motor	250W 3:1 Gear Motor	400W 3:1 Gear Motor	400W 3:1 Gear Motor
Speed	60-600RPM	60-600RPM	60-600RPM
Stirrer	Stainless steel with PTFE Coating 12mm	Stainless steel with PTFE Coating 16mm	Stainless steel with PTFE Coating 25mm
Maximum Temperature	200℃	200℃	200℃
Minimum Temperature	-120℃	-120℃	-120℃
Maximum Vacuum	-0.098Mpa	-0.098Mpa	-0.098Mpa
Voltage	220V	220V	220V
Package	Craft Box	Craft Box	Craft Box

Note: Jacketed glass reactor can be for heating and cooling ,but you need purchase heater or cooler additionally



科瑞仪器

>>> Kori instrument

Xinyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xinyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Turnkey setup System: Optional 1

Jacketed glass reactor with Heat circulating bath ---- Suitable for heating and can heat up to 200°C (Customized 300°C)



科瑞仪器

>>> Kori instrument





科瑞仪器

>>> Kori instrument

Xingyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xingyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Turnkey setup System: Optional 2

Jacketed glass reactor with recirculating chiller -----Suitable for cooling and Minimum temperature as you need



科瑞仪器

>>> Kori instrument





科瑞仪器
>>> Kori instrument

Xingyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xingyang, Zhengzhou City, 450100, China

Tel: 0086-371-86573886

Fax: 0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Turnkey setup System: Optional 3

Jacketed glass reactor with recirculating chiller -----Suitable for heating and cooling

Maximum temperature : 200°C

Minimum temperature: as you need



科瑞仪器
>>> Kori instrument





科瑞仪器
Kori instrument

Xingyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xingyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Full Setup Turnkey system : Vacuum pump +Jacketed glass reactor+chiller+heating cooling machine





科瑞仪器

>>> Kori instrument

Xinyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xinyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Jacketed glass reactor with Ex-proof





科瑞仪器
Kori instrument

Xinyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xinyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Technical data **10L 20L 30L 50L Glass Reactor (Ex-proof)**

Model	SF-10LEX	SF-20LEX	SF-30LEX	SF-50LEX
Size	10L	20L	30L	50L
Motor	180W Ex-proof motor	180W Ex-proof motor	180W Ex-proof motor	180W Ex-proof motor
Speed	60-600RPM	60-600RPM	60-600RPM	60-600RPM
Stirrer	Stainless steel with PTFE Coating 12mm	Stainless steel with PTFE Coating12mm	Stainless steel with PTFE Coating 12mm	Stainless steel with PTFE Coating 12mm
Maximum Temperature	200℃	200℃	200℃	200℃
Minimum Temperature	-120℃	-120℃	-120℃	-120℃
Maximum Vacuum	-0.098Mpa	-0.098Mpa	-0.098Mpa	-0.098Mpa
Voltage	220V	220V	220V	220V
Package	Craft Box	Craft Box	Craft Box	Craft Box

Note: Jacketed glass reactor can be for heating and cooling ,but you need purchase heater or cooler additionally



科瑞仪器

>>> Kori instrument

Xinyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xinyang, Zhengzhou City, 450100, China

Tel: 0086-371-86573886

Fax: 0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130





科瑞仪器
Kori instrument

Xingyang Kori Instrument Factory (Henan Zhongyi Kori Equipment Co.,Ltd)

Address: Mugou Village, Gaoshan District, Xingyang,Zhengzhou City,450100,China

Tel:0086-371-86573886

Fax:0086-371-64877188

www.korichina.com

E-mail: linda@korichina.com/info@korichina.com

Mobile/Wechat/Whatsapp : 0086 13683828130

Technical data 100L 150L 200L **Glass Reactor (Ex-proof)**

Model	SF-100LEX	SF-150LEX	SF-200LEX
Size	100L	150L	200L
Motor	370W Ex-proof Motor	750W Ex-proof Motor	750W Ex-proof Motor
Speed	60-600RPM	60-600RPM	60-600RPM
Stirrer	Stainless steel with PTFE Coating 12mm	Stainless steel with PTFE Coating 16mm	Stainless steel with PTFE Coating 25mm
Maximum Temperature	200℃	200℃	200℃
Minimum Temperature	-120℃	-120℃	-120℃
Maximum Vacuum	-0.098Mpa	-0.098Mpa	-0.098Mpa
Voltage	220V	220V	220V
Package	Craft Box	Craft Box	Craft Box

Note: Jacketed glass reactor can be for heating and cooling ,but you need purchase heater or cooler additionally

Welcome you to visit our factory ;)