



## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel:+86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### Choicy Beauty Equipment Series

### -- Portable Cryolipolysis Body Shaping and Slimming System



**Model: CB-702C**

### Product Description





## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel:+86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### System Characteristics

1. 8 inch colorized touch screen control, simulation modeling of dynamic 3D interface, human-computer interaction design unique, simple and convenient operation
2. convenient operation
3. Patent design of separation technology, can effectively prevent the skin in a long time below 5 degrees Celsius.
4. Exclusive development of freezing point reduction and surface active dosage form antifreeze film, can effectively prevent skin hurt;
5. Automatic water adding and monitoring function, can be filled with water tank and hand channel at one time
6. Automatic cleaning and monitoring function to well protect the handpiece
7. The intelligent monitoring system can ensure the normal operation of the whole system
8. Intelligent alarm system can quickly and effectively respond to the abnormal situation
9. Super intelligent self inspection system, can quickly and effectively find out the problem part in the abnormal situation

### Treatment Range

1. Fat reduce, loss weight, Body Slimming
2. Tighten skin, anti-aging
3. Improve blood circulation and metabolism
4. Smooth scar and wrinkles
5. Smooth striae gravidarum
6. Add skin moisture
7. Enhance fibroblast elastic tissue

## Before and After Treatment



## Technical Specifications

Input voltage	AC110V/220V 50-60Hz
Wavelength	Red (630nm) Blue (470nm)
Screen	8 inch wide color touch screen
Cryo	+5°C to -12°C
Vacuum pressure	10Kpa-80KPa
Cooling liquid	Purified water or special coolant
Ambient temperature	5°C to +40°C
Relative humidity	≤80%
Atmospheric pressure	90KPa-110KPa
Packing size	67X47X78 cm

## NOTE:

The detailed instruction will be sent with the equipment.