



## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel:+86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### Choicy Beauty Equipment Series -- Phototherapy System



**Model: CB-903A**



## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province,  
China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel:+86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### System Characteristics

1. LED SYSTEM is a phototherapy machine which professional focuses on increase skin tissue synthesis of collagen and elastin fibres, to produce healthier, faster healing, younger looking and glowing skin.
2. LED technology that generates photons at different frequencies to penetrate into deeper layers of tissue. Different wavelengths (colors blue, infrared, red) of light will have different effects on the skin.
3. LED SYSTEM with Multi functions which is recommended for improving oily skin, removing marks and excess pigmentation, diminishing wrinkles, skin flaccidity, expression lines and orange-peel skin, for firming skin, pain relief.
4. Built-in intelligent treatment can solve problems in a targeted manner
5. Electrically controlled lifting rod height, 270 degree rotating treatment lamp board, easier focus on treatment area

### Treatment Range

1. Remove wrinkles, and pigments on facial skin
2. Skin Rejuvenation and anti-aging
3. Recovery of injured skin:
4. Smooth Lymph and Eliminate Edema.
5. Improve oily skin, to reduce the chances of getting acne.
6. Pain relief

## Before and After Treatment



## Technical Specifications

Model	CB-903A
LED Dominant Wavelength	RED: 650nm IR: 810nm BLUE: 450nm
Laser type	3 colors LED
Treatment panel	4 pieces
Energy	110mw/cm2
Cooling Method	Strong wind Cooling
Screen	8 Inch LCD Color Touch Screen
Input Power	800W
Power Supply	AC110V/8A/60Hz AC220V/4A/50Hz
Working time	1-60min(adjustable)

## NOTE:

The detailed instruction will be sent with the equipment.