

Add: Room 1203 No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

Choicy Beauty Equipment Series

-- Hydra-Revival Skin Care



Model: CB-1101A



Add: Room 1203 No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

System Characteristics

- 1. Hydra -revival skin care machine is the "skin cleansing skin care" function combined with a comprehensive beauty machine.
- 2. A multifunctional machine that uses a variety of sterile small bubbles, revival-oxygen, RF and ultrasound therapy, four kinds of techniques to treat a variety of skin problems, including acne, blackheads, skin aging, dehydration, sensitive skin and fine lines.
- 3. Noninvasive non-surgical: non-invasive skin tightening, no hurt, safe and effective.
- 4. Suitable to wider population, all types of skin surgery can be used safely before surgery.

Treatment Range

- 1. Aqua cleaning handle-piece deep cleansing for skin and removing the mites.
- 2. Revival & oxygen handle-piece— deep cleansing, eliminate the aging of keratinocytes, remove the sebum, grease residue, as the same time put the skin into the oxygen and remove the yellow skin.
- 3. RF handle-piece—to improve fine lines, and skin tightening all round.
- 4. Ultrasound import handle-piece— put the skin into the nutrients, and provide lasting nutrition, to make the skin moist and delicate shiny.

Technical Specifications

Model	CB-1101A
-------	----------



Add: Room 1203 No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

Handpiece quantity	4 pcs	
Handpiece type	Aqua cleaning handpiece, Revival & oxygen handpiece, RF handpiece,	
	Ultrasound import handpiece	
Handle connection	Quick connection	
Aqua cleaning handpiece		
Procrastination	Super	
Size	Big (19mm*15mm) small (13mm*905mm)	
Vacuum pressure	-30Kpa~-70Kpa	
Revival & oxygen handpiece		
Tips	Vibration massage tip, Oxygen powered tip (Optional)	
Diameter	28mm	
RF handpiece		
RF Energy	10—100	
Frequency	1MHz	
Max Power	10W	
Peak Energy	200+20V	
Ultrasound import handpiece		
Ultrasound	10-100	
Frequency	1MHz	
Max power	10W	
Screen	10.4 inch touch screen	
Input power	100W (Max)	
Output power	90W (Max)	

Before and After Treatment



Add: Room 1203 No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province China

E-mail: info@choicybeauty.com

Tel:+86-20-38263199

Website: www.choicybeauty.com









NOTE:

The detailed instruction will be sent with the equipment.