

Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

Choicy Beauty Equipment Series

-- Dr. Pen System



Model: A6

Product Description



Optional Needles




Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

Needle



Introduction

Minimised epidermal damage

Less pain and erythema as high speed with vertical direction

More easy to operate at the small or curved areas

Adjust micro needle length (0.25~3.0mm)

Prevent cross infection, adopt with disposable needles

More quickly and efficient than manual Micro-needle

System Characteristics

- 1. Atomatic and charging
- 2. Adjustable Needle length from 0.25mm to 3.0mm
- 3. Adjustable 6 level speed control from 1800rpm to 7000rpm
- 4. 12 / 9 needles Tip using 30 gauge Acupuncture Needles to reduce discomfort
- 5. Digital Motor in Hand-piece
- 6. Highest quality for medical professionals(Aluminum Housing)
- 7. Optimal absorption of products
- 8. Compact and ergonomic design

Advantages

Choicy Beauty

Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

- A. Visual and noticeable improvement in a few days from the first treatment.
- B. Short healing time.
- C. No permanent damage or injury.
- D. No bruising, infection, discoloration or other complications.
- E. Extremely high absorption of any active ingredients.
- F. Almost all skin types can be treated.
- G. thickens the dermis without damaging or removing epidermis.

Treatment Range

- 1. Anti-Puffiness.
- 2. Scar removal including acne scar removal or treatment.
- 3. Acne Scar Removal, Acne Scar Treatment.
- 4. Anti Wrinkle Treatment, Anti Aging.
- 5. Hair Loss Treatment or Hair Restoration.
- 6. Surgical Scar Removal.
- 7. Hyper pigmentation treatment.

Specification:

Туре	derma Pen
Power supply	charging or battery
Adapter	4.5v-500mA
Vibrating frequency	1800rpm - 7000rpm
Needle length	0.25~3.0mm
Needle material	nano wafer /imported medical stainless steel
Needle number	12/9/36/42 needles, nano pins, silicon pins for facial treatment 1/3/5/7 needles for tattoo
Body material	aluminum housing



Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

Body color	Silver/red/blue
Size	130mm
Weight	65g

NOTE:

The detailed instruction will be sent with the equipment.