



Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou,
Guangdong Province, China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com


Choicy Beauty Equipment Series

-- RF Fractional laser & non-invasive resurfacing system



Model: CB-807B

Product Description



System Characteristics

1. More tiny and even-energy dot

Minimum 120 μ m dot makes treatment goes into dermis and especially reducing the risk of producing pigmentation after treatment, where it can produce significant dermal impact.

2. Fast treatment

Big treatment area and high pulse energy, can finish the whole face treatment within only 10 minutes

3. Short downtime

Patients can quickly resume their normal activities.

4. Safe and effective for all skin types

5. Air cooling without water cooling, avoids more troubles causing by complicated water cooling.

Decrease the failure rate.

Treatment Range

1. wrinkle removal

2. Skin resurfacing

3. Acne scar removal

4. Scald and surgical scar recovery

- 5. Pigmentation removal
- 6. Photo - damage - skin recovery
- 7. Melanotic nevus removal

Before and After Treatment



Technical Specifications

Laser type	RF laser tube imported from USA
Laser mode	Super pulse RF laser
Laser wavelength	10600nm ± 5nm
Beam-divergence angle	4.5 mard
Laser output power	1-30w (adjustable), 1w/step, deviation ± 10%
Indicated light wavelength	650 ± 5nm(diode)
Cooling system	Strong wind cooling
Screen	10.4 inch color touch screen
Scan range	20mm*20mm
Input voltage	AC110V/15A/60Hz, AC220V/8A/50Hz
Net weight	55kg
Dual mode	
Normal mode	
Output type	Single pulse, multiple pulse, continuous pulse, magical pulse
Pulse width	10ms-500ms(adjustable)



Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou,
Guangdong Province, China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

Pulse interval	0.2-1.5s(adjustable)
Scan mode	
Output type	Graphics scanning
Optional graphics	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Laser pulse width	1-50ms (adjustable)
Power	1-30w (adjustable)

NOTE:

The detailed instruction will be sent with the equipment.