

## **Choicy Beauty Equipment Series**

### **-- Nano Laser**



**Model: CB-306A**

### **System Characteristics**

1. Pigment removal speed is faster
2. Pigment clear rate is higher
3. Don't damage the normal tissue
4. Refused to stimulate melanin active again risk

### **Treatment Range**

1. Refractory chloasma, coffee spots, age spots, freckles, sunburn, nevus of Ota, etc
2. Skin whitening, wrinkles removal
3. Acne scars

## 4. All colors of tattoo removal

**Technical Specifications**

Laser wavelength	1064nm/532nm
Laser type	Solid laser
Energy	1064nm: 300-1200mJ (adjustable) 532nm: 160-600mJ (adjustable)
Pulse width	6ns ± 2ns
Beam-divergence angle	1 mrad
Pulse frequency	1-10Hz (adjustable)
Spot size	1-10mm (adjustable)
Cooling system	Enclosed internal circulation water cooling, external circulation of strong wind cooling
Indicated light wavelength	650nm
Indicated light intensity	1-5 (adjustable)
Screen	10 inch touch screen
Input power	2200w (AC220V / 10A, AC110V/20A)
Net weight	50kg

**Before and After Treatment**

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province,  
China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel:+86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

Traumatic tattoo removal in just  
one treatment



After two treatments



After four treatments



After five treatments



## NOTE:

The detailed instruction will be sent with the equipment.