

# Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1502 Hongmian Business Center, No.10 Hualei Road, Liwan District, Guangzhou, Guangdong Province, China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

### **Choicy Beauty Equipment Series**

## -- E-Light (IPL+RF) Skin Rejuvenation System



Model: CB-601A

#### **Product Description**







### Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1502 Hongmian Business Center, No.10 Hualei Road, Liwan District, Guangzhou, Guangdong Province, China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

#### **System Characteristics**

- 1. Treatment result is obviously enhanced on clinic, and the various skin problems like wrinkle and telangiectasia are greatly improved.
- 2. Multi-features of E-light machine lead to the machine as the best choice for patient: Comfortable treatment, better effect, faster results and minimal side effects, etc.
- 3. Adapting to the wider population, it can be used for the treatment of skin Type IV and skin type V.
- 4. Easy operation, you only need to adjust the energy of IPL and RF individually.
- 5. Safe to all kinds of skin types.
- 6. Pre-cooling to the epidermis, this can result in temperature difference between epidermis and dermis, so as to achieve the best RF penetration effect;
- 7. Pre-cooling the skin, the energy lunched by RF act on the pre-warmed aimed-tissue (such as veins, hair follicles, melanin, etc.)
- 8. The treatment head of E-light mainly combined with three parts: Bipolar radio frequency + low-energy pulse light energy + epidermal cooling.

### **Treatment Range**

- 1. Hair removal: suitable for both dark and light color skin, especially for light color skin.
- 2. Skin rejuvenation: remove wrinkles, speckles, age spot, sun-induced freckles
- 3. Skin tightening, skin lifting, improves flexibility and tone of skin.
- 4. Acne removal: improve the situation of oily skin.
- 5. Vascular removal, rosacea treatment, red veins removal



# Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1502 Hongmian Business Center, No.10 Hualei Road, Liwan District, Guangzhou, Guangdong Province, China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

#### **Before and After Treatment**



### **Technical Specifications**

Model	CB-601A
RF power	5-100W
RF output time	0.2-3.5s
Spot size	8x34mm
Standard handpiece(wavelength)	560-1200nm(skin rejuvenation),695nm-1200nm(hair removal)
Optional handpiece(wavelength)	400-1200nm,530-1200nm,595-1200nm
Energy density	10 - 60J/cm2
Pulse number	1 - 6pcs
Pulse delay	5 - 100ms
Pulse width	1 - 10ms
Operate Interface	8 Inch Color Touch Screen
Cooling system	Semiconductor & Water & air cooling system
Power Supply	AC110V/15A/60Hz AC220V/8A/50Hz
Frequency	1-3Hz
Dimensions	52cm*52cm*95cm
Weight	55KG

### **NOTE:**

The detailed instruction will be sent with the equipment.