

Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203 No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

Choicy Beauty Equipment Series

-- Vacuum Body Slimming Physiotherapy



Model: CB-502A

Product Description





Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203 No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

System Characteristics

This machine adopts vacuum technology, the early application of vacuum in medicine, to remove scar and recovery after liposuction. With the progress of technology, vacuum is applied to the deep massage, promoting cells do aerobic exercise, with extraordinary detoxification health care efficacy and also has good effect on improving stretch marks tissue, focus on the "health body shaping"

- 1. Vacuum and roller system, equipped with 12 modes
- 2. RF: high-power RF output system
- 3. Near-infrared laser wavelength system: semiconductor laser technology at 650nm
- 4. New frequency conversion air pressure control technology, constant air pressure
- 5. Rapid relief mode around the face and eyes

Treatment Range

- 1. Postpartum repairing and body shaping
- 2. Improving stretch marks, remodeling the body and the waistline becomes smaller after liposuction
- 3. Systemic or partial obesity, wrinkle removal, fat-soluble and skin tightening and body shaping (arms, legs and shoulder)
- 4. Relieving muscle pain, systemic physical therapy, etc
- 5. Pregnancy maintenance, improve stretch marks to varying degrees
- 6. Scar repair, can effectively alleviate scar tissue
- 7. Rapid treatment around the eyes, relieve dark circles and repair eye bags
- 8. Body massage, can massage with muscle texture



Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203 No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province China

E-mail: info@choicybeauty.com Tel:+86-20-38263199 Website: www.choicybeauty.com

Technical Specifications

Operation mode	Pulse mode
Pulse width	0.1s - 9.9s
Vacuum pressure	80kpa-10kpa (60.8Hg-7.6Hg) (Absolute value)
	20kPa - 90kPa (15.2Hg- 68.4Hg) (Relative value)
Roller speed	12-33rpm
Roller working mode	4 types
RF Frequency	1MHz
RF energy density	60J/cm3 (max)
Near-infrared laser wavelength	650nm
Near-infrared laser power	10W(max)
Hand-piece quantity	4PCS
Rated voltage	AC 230V(+/-10%) 50Hz,AC 110V(+/-10%) 60Hz
Screen of the machine	TFT color touch screen
Screen of vacuum Hand-pieces	Big hand-piece: 2.4"
	Small hand-piece:1.9"
Real-time Online Safety Detection	

Before and After Treatment



NOTE:

The detailed instruction will be sent with the equipment.