



## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel:+86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### Choicy Beauty Equipment Series

### --Fractional laser & non-invasive resurfacing system



**Model: CB-803A**

#### Product Description





## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel:+86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### System Characteristics

1. 3D fractional laser technology makes treatment be much safer and more effective.
2. High-power and super-stable laser offers higher energy.
3. Multi-fractional modes
4. Adjustable rectangular pulse width control technology precisely controls pulse and strictly controls heat transfer.
5. Personalized touch screen suits to different customer requirement
6. Non-invasive treatment,short down time and comfortable treatment.

### Treatment Range

#### 1. Skin Rejuvenation

The skin is micro exfoliated under fractional mode, and heat stimulation of fractional laser makes collagen renew to get more tender and more elastic skin.

#### 2. Wrinkle Removal

Remove wrinkles around mouth, eyes, forehead etc

#### 3. Scar Treatment

Remove surgical scars and acne scars

#### 4. Freckle Removal

Dilute various freckles like chloasma, freckles, coffee spots, age pigment and the actinic keratoses.

#### 5. Remove skin neoplasms, Syringoma, sebaceous, gland tumor, Becker nevus, and seborrheic keratosis.

#### 6. Improve skin texture, shrink pores, make the skin much smoother and tenderer

#### 7. Tightening face and neck, with the function of face-lift and enhancement

#### 8. Whitening face and remove pigmentation



# Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: info@choicybeauty.com    Tel:+86-20-38263199    Website: www.choicybeauty.com

## Clinical Applications

1. To improve skin texture, shrink pores, make the skin smooth and tender.
2. To tighten Face and neck, with the function of face-lift and enhancement.
3. To remove fine wrinkles around eyes and mouth.
4. To treat surgical and acne scars.
5. To whiten face and remove pigmentation.

## Technical Specifications

Laser type	Sealed off CO <sub>2</sub> laser
Laser wavelength	10.6 μ m
Laser mode	Low order mode
Output power	1-30W (adjustable)
Focus spot diameter	0.4mm
Lens focal distance	F=50mm/F=100mm
Delivery system	7-joint Articulated arm
Working modes	Continuous pulse, single pulse, repeat pulse and super pulse
Pulse duration	1-500ms
Cooling system	Closed Loop Circulating Water
Power supply	110V/8A 60Hz 220V/4A 50Hz
Screen	8 inch color touch screen
Scan range	20*20mm
Scan graphics	rectangle, ellipse, triangle, and parallelogram
Scan density	0.2-2mm/dot
Gross weight	45kg



## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel:+86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### Before and After Treatment



#### NOTE:

The detailed instruction will be sent with the equipment.