



**Guangzhou Huaxi C&T Co.,Ltd.**

Add: Room 1502 Hongmian Business Center, No.10 Hualei Road, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel: +86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

## **Choicy Beauty Equipment Series**

### **-- RF Vacuum Slimming System**



**Model: CB-501A**

### **Product Description**





## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1502 Hongmian Business Center, No.10 Hualei Road, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel: +86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### System Characteristics

#### 1. Non-invasive body slimming

Optimum touch through multi-angle operation: Using electrode vertically & crossly

#### 2. Mechanical operation

Massage: Strengthen the RF & heating penetration depth, effectively arrived the fatty layer.

Vacuum: Closer contacting the treatment target, increasing blood circulation, and promoting cell metabolism.

#### 3. Thermal effect

Bipolar RF- Intensively heating fat cells layer, effectively removing stretch marks.

IR-Adjusting skin elasticity, guiding the RF energy to the connective tissue deeply.

#### 4. Multi-mode hand-piece

Big hand-piece: Treat in big area, such as belly, waist, thighs, hips etc.

Small hand-piece: Treat in small area, such as arms, neck, belly, face and thin fat area etc.

### Treatment Range

1. Shape the body.
2. Cellulite reduction.
3. Body shaping by reducing cellulite.
4. Reduce measure after fat dissolving.
5. Pretreatment by measure-reduction
6. Eliminate the stretch marks
7. Skin tightening
8. Massage



## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1502 Hongmian Business Center, No.10 Hualei Road, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel: +86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### Technical Specifications

Operation mode	Pulse mode
Vacuum pressure	80kpa-10kpa (Absolute value)
Roller speed	12-33rpm
Roller working mode	4 types 1. Roll in mode(two wheels of the big hand-piece roll toward inside) 2. Roll out mode(two wheels of the big hand-piece roll toward outside) 3. Roll up mode(two wheels of the big hand-piece roll up or down in different speed) 4. Normal mode (two wheels of the big hand-piece roll up or down in same speed)
RF Frequency	1MHz
RF energy density	60J/cm <sup>3</sup> (max)
Near-infrared laser wavelength	650nm
Near-infrared laser power	20W(max)
Hand-piece quantity	4PCS
Rated voltage	AC 220V(+/-10%) 50Hz, AC 110V(+/-10%) 60Hz
Rated power	1500W
Rated current	6.8A
Working ambient temperature	0-30℃
Storage temperature	0-45℃
Relative humidity	≤ 80%
Screen of the machine	10 inch TFT color touch screen
Screen of vacuum Hand-pieces	Big hand-piece: 2.4" Small hand-piece: 1.9"



## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1502 Hongmian Business Center, No.10 Hualei Road, Liwan District, Guangzhou,  
Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel: +86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### Before and After Treatment



### NOTE:

The detailed instruction will be sent with the equipment.