



ENDO LINEAR CUTTER STAPLER AND RELOADING UNIT

www.cnweipumedical.com
Email: lisa@cnweipumedical.com
Tel: + 86 198 5020 8586

Direction for Use

The product is applied to resection, transection and anastomosis of lung, bronchial tissue, stomach and intestine in open or endoscopic surgery.



Features

-The stapler handles with an ergonomic design (with rubber grip), can be one-handed grasping which provides a comfortable and easy control;

-Slimmer fixed anvil allows for enhanced clamping force, improves maneuverability of the stapler around target tissue;

-Varied height staples design can reduce tissue compression stress on outer staple;

-Gentle tissue compression before firing minimizes tissue trauma and supports re-positioning;

-Cartridge lock out safety mechanism prevents use of fired reloading units;

-One universal handle can fitted with all reloading units in different size.



Stapling System




Specification Details	
Specification	Shaft Length (MM)
WPYQB-60	60
WPYQB-160	160
WPYQB-260	260

Reloading Unit



Specification	Staple Height(MM)
WPJQ45M-D	3 / 3.5 / 4
WPJQ45H-D	4 / 4.5 / 5
WPJQ60L-D	2 / 2.5 / 3
WPJQ60M-D	3 / 3.5 / 4
WPJQ60H-D	4 / 4.5 / 5
WPJQ4525-D	2.5
WPJQ6035-D	3.5

WPJQ45L-W	2 / 2.5 / 3
WPJQ45M-W	3 / 3.5 / 4
WPJQ60L-W	2 / 2.5 / 3
WPJQ60M-W	3 / 3.5 / 4
WPJQ4525-W	2.5
WPJQ6035-W	3.5



Application Reference

Reload Specification	Staple Height	Tissue Thickness (MM)
WPJQ45L-D (Brown)	2/2.5/3	0.88~1.80
WPJQ45M-D (Purple)	3/3.5/4	1.50~2.25
WPJQ45H-D (Black)	4/4.5/5	2.25~3.00
WPJQ60L-D (Brown)	2/2.5/3	0.88~1.80
WPJQ60M-D (Purple)	3/3.5/4	1.50~2.25
WPJQ60H-D (Black)	4/4.5/5	2.25~3.00
WPJQ4525-D (White)	2.5	0.88~1.80
WPJQ6035-D (Blue)	3.5	1.50~2.25
WPJQ45L-W (Brown)	2/2.5/3	0.88~1.80
WPJQ45M-W (Purple)	3/3.5/4	1.50~2.25
WPJQ60L-W (Brown)	2/2.5/3	0.88~1.80
WPJQ60M-W (Purple)	3/3.5/4	1.50~2.25
WPJQ4525-W (White)	2.5	0.88~1.80
WPJQ6035-W (Blue)	3.5	1.50~2.25

