

Choicy Beauty Equipment Series

-- Fractional CO2 Vaginal Tightening Laser



Model: CB-812A

Product Description



System Characteristics


1. Tightens the vagina and rejuvenates natural flexibility
2. Vaginal remodeling pain free procedure with long-lasting result
3. Greater sexual pleasure and desire.
4. Increase quantity and quality of orgasms achieved during intimacy.
5. Almost non-painful, no need for anesthetics during treatment
6. Improves vaginal dryness and stress urinary incontinence
7. Very short treatment time, just about 15-20min at clinic

Treatment Range

1. Tighten vagina
2. Clean the vagina
3. Moist vagina
4. Improve skin texture & skin tone
5. Skin resurfacing
6. Collagen contraction and remodeling
7. Pore reduction
8. Pigmentation & age spots removal
9. Photo aged skin and fine lines treatment
10. Scars removal – chicken pox, post acne, post traumatic scars

Before and After Treatment


Vaginal Canal Tighten



Vaginal Tighten



Technical Specifications

Laser type	Sealed off CO ₂ laser
Laser wavelength	10.6 μ m
Laser mode	Low order mode
Output power	1-30W (adjustable)
Focus spot diameter	0.4mm
Lens focal distance	F=50mm/F=100mm
Delivery system	7-joint Articulated arm
Working modes	Continuous pulse, single pulse, repeat pulse and super pulse
Pulse duration	1-500ms
Cooling system	Closed Loop Circulating Water
Power supply	110V/8A 60Hz 220V/4A 50Hz
Screen	8 inch color touch screen
Scan range	20*20mm
Scan graphics	
Scan density	0.2-2mm/dot
Gross weight	45kg

NOTE:

The detailed instruction will be sent with the equipment.