

**Choicy Beauty Equipment Series**  
**-- Youthful RF System**



**Model: CB-1303A**

**System Characteristics**

1. Monopolar RF technology, deep heating
2. Precise skin temperature detection
3. Quickly plug and unplug, intelligent recognition of the head
4. 10mm to 30mm has different treatment head, for all parts

## Technical Characteristics

1. It can be used fixed-point treatment modes, each part reaches the preset temperature and stop
2. Accurate temperature feedback, real-time monitoring of epidermis temperature
3. Intelligent identification of the plug-in treatment head, automatic ally the output parameters

## Treatment Range

1. Skin tightening and improve flabby
2. Cellulite fat treatment
3. Promote metabolism and blood circulation
4. Improve stretch marks
5. Muscle release physiotherapy

## Before and After Treatment



Before

5 treatments

3 Month follow-up



Before

6 treatments

3 Month follow-up



## Guangzhou Huaxi C&T Co.,Ltd.

Add: Room 1203, No.66 Huadi South Ave, Liwan District, Guangzhou, Guangdong Province, China

E-mail: [info@choicybeauty.com](mailto:info@choicybeauty.com) Tel:+86-20-38263199 Website: [www.choicybeauty.com](http://www.choicybeauty.com)

### Technical Specifications

RF frequency	1MHz
Maximum output energy	200w
Treatment handpiece area	Body treatment handpiece: 30mm Neck treatment handpiece: 25mm Facial treatment handpiece: 18mm Eyes treatment handpiece: 10mm
Display	10 inch color screen touch screen
Power supply	220+-10%, 50/60 HZ 5A
Weight	20kg

### NOTE:

The detailed instruction will be sent with the equipment.